Medtronic
MINIMED

SOF-SENSOR™

**GLUCOSE SENSOR • CAPTEUR DE GLUCOSE
GLUKOSESENSOR • SENSOR DE GLUCOSA
GLUCOSESENSOR • SENSORE DELLA GLICEMIA
GLUKOSESENSOR • GLUKOOSISENSORI
GLUKOSENSOR • ΑΙΣΘΗΤΗΡΑΣ ΓΛΥΚΟΖΗΣ
GLUKOSESENSOR • SENSOR DE GLICOSE
CZUJNIK GLUKOZY • SENZOR PRO MĚŘENÍ
GLUKÓZY • GLÜKÓZÉRZÉKELŐ
GLUKÓZOVÝ SENZOR • SENZOR ZA GLUKOZO
GLIKOZ SENSÖRÜ • СЕНСОР ГЛЮКОЗЫ
جهاز استشعار الغلوكوز • חישן סוכר**

R_x Only

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






Indications For Use

The Glucose Sensor is intended for use with Medtronic MiniMed glucose sensing systems to continuously monitor glucose levels in persons with diabetes.

Warnings

Failure to follow instructions may result in pain or injury.

The Sensor may create special needs regarding your medical conditions or medications. Please discuss these conditions and medications with your doctor before using the Sensor.

Bleeding, swelling, irritation and/or infection at the insertion site are possible risks associated with inserting the Sensor and may result from improper insertion and maintenance of insertion site.

After insertion, wait 5 minutes before attaching the transmitter to watch for bleeding at the insertion site. If bleeding occurs, apply steady pressure using a sterile gauze or clean cloth for up to 3 minutes.

If bleeding stops:

1. Attach the sensor to the transmitter.

If bleeding does NOT stop:

CAUTION: DO NOT attach the sensor to the transmitter.

1. Remove the sensor and discard.
2. Check the site for redness, bleeding, irritation, pain, tenderness or inflammation and treat accordingly.
3. Insert a new sensor in a different location

The Sensor is sterile and non-pyrogenic, unless the package has been opened or damaged. Do not use any Glucose Sensor if its sterile package has been previously opened or damaged. Always inspect packaging for damage prior to use.

If Sensor is not securely placed in Sen-serter[®] prior to insertion, pain or minor injury may occur.

Remove needle guard before inserting Sensor. If you develop inflammation, redness, soreness or tenderness at insertion site, or if you experience unexplained fevers, remove Sensor. Check site often for these conditions and to ensure Sensor is still in place, especially before going to bed and upon waking.

Precautions

Always wash hands with soap and water before opening the Sensor package and handling Sensor. After opening the package, avoid touching any Sensor surfaces that will come in contact with the body, i.e., Sensor, needle, adhesive surfaces and tape.

Before inserting Sensor, clean insertion site with a topical antimicrobial solution, such as isopropyl alcohol. Allow area to dry.

Do not insert Sensor through tape, and do not use sticky preps before insertion.

Establish a rotation schedule for choosing each new Sensor site. Choose a site with adequate subcutaneous fat for Sensor insertion. The lower abdomen and upper, outer quadrant of the buttocks are most frequently used. Avoid sites that are constrained by clothing, scar tissue,

accessories or subjected to rigorous movement during exercise. Also, avoid sites around belt or waistline or 2 inches (5.08 cm) around navel.

Make sure Sensor insertion site is at least 3 inches (7.62 cm) from insulin pump infusion site or manual injection site. When replacing Sensor, select new insertion site at least 2 inches (5.08 cm) from previous site.

After inserting Sensor and attaching it to the transmitter, cover the Sensor site with sterile/occlusive dressing, such as IV3000[®] (MMT-174).

Dispose of used Sensors and introducer needles in a sharps container after single use. Do not clean or resterilize, and do not recap the needle with Sensor needle cover.

Reagents

The Sensor contains two biological reagents: Glucose Oxidase and Bovine Serum Albumin (BSA). Glucose Oxidase is derived from *Aspergillus niger*, purified and dried according to Type VII-S guidelines and cross linked to BSA, purified and dried albumin fraction V, derived from pasteurized bovine serum with gluteraldehyde. Since less than 0.4 µg of glucose oxidase and less than 0.7 µg of BSA is used to manufacture each Sensor, the risks of tissue reactions and viral transmission are considered minimal.

Storage and Handling

Sensors can be stored at room temperature between +36°F to +80°F (+2°C to +27°C). For temperatures greater than +80°F (+27°C), sensors will require cooled storage at temperatures not lower than +36°F (+2°C). Prior to using a sensor stored in a cool environment, allow the sensor

package to warm up to room temperature for about 15 minutes (at 5% to 95% relative humidity) before opening of the sensor package and use, as to prevent condensation.

CAUTION: Do not freeze the sensor. Sensors should not be stored in direct sunlight, extreme temperatures, or humidity.

NOTE: Discard Sensor after “Use Before” date on label, or if the package is damaged or the seal is broken.

See your Sen-Serter User Guide for detailed information concerning your Sen-Serter.

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- EN** The Sensor and Sen-serter are covered by the following U.S. Patents: 5,390,671; 5,391,250; 5,568,806; 5,586,553; 5,777,060; 5,779,665; 5,786,439; 5,851,197; 5,882,494; 5,954,643; 6,093,172; 6,293,925; 6,462,162; 6,520,326; 6,607,509; DES 433,755. The Glucose Sensor and Sen-serter are covered by other/pending patents. ©2004, Medtronic Minimed, all rights reserved. CGMS[®], MiniMed[®] and Sen-serter[®] are registered trademarks of Medtronic MiniMed. Guardian[™] is a trademark of Medtronic MiniMed. IV3000[®] is a registered trademark of Smith&Nephew, Inc.
Medtronic B.V.