

Please use the following checklists to make sure you have what you need during your in-person training sessions.

## For your training sessions

Always bring your Getting Started Guides to any in-person training sessions.

- Getting Started Guide Packet
  - GETTING STARTED WITH YOUR MINIMED® 670G INSULIN PUMP
  - GETTING STARTED WITH MINIMED® 670G CONTINUOUS GLUCOSE MONITORING
  - GETTING STARTED WITH YOUR MINIMED® 670G AUTO MODE
  - Your letter and this checklist from the Training Packet

## Insulin pump training

### If you are new to pump therapy

Make sure you stop using your long acting insulin the day before training. Talk to your healthcare professional to determine treatment options during this time.

- MiniMed 670G insulin pump
- The blue and black Medtronic bag with all contents:
  - Reservoirs
  - Infusion sets
  - Serter
  - CONTOUR® NEXT LINK 2.4 meter
- Three I.V. PREP™ Antiseptic Wipes (or alcohol wipes)
- A vial of rapid-acting insulin U-100

The reservoir for your insulin pump cannot be filled with insulin from an insulin pen. You must have a vial of insulin as prescribed by your healthcare professional.

- Pump Initiation Settings/Prescriber's Instructions to Patient form completed and signed by your healthcare professional with your initial pump settings (your certified product trainer may already have this)
- Other
  - Glucose tablets or some form of fast-acting carbohydrate
  - A snack

## For CGM training

CGM training may be done alone, or with Auto Mode. Either way, you will need to bring the following with you for your CGM training.

- ❑ MiniMed 670G insulin pump
- ❑ Guardian® Link transmitter kit:
  - Fully charged Guardian Link transmitter
  - Gray charger (with AAA battery)
  - One-press Serter device
- ❑ Two Guardian glucose sensors
- ❑ Enlite® overtape located in a small plastic bag in the Guardian sensor box

If you bring a box of sensor to the training, be sure to not leave the sensors in a hot car. Heat can cause damage to the sensors.

- ❑ Alcohol wipes
- ❑ CONTOUR® NEXT LINK 2.4 meter
- ❑ Test strips
- ❑ Lancets
- ❑ Other
  - Glucose tablets or some form of fast-acting carbohydrate
  - A snack

## For Auto Mode training:

If you have already started using your MiniMed 670G insulin pump with CGM, do the following:

- Calibrate your sensor prior to training
- Bring your CONTOUR NEXT LINK 2.4 meter
- Ask your trainer if you should bring a snack