



# TRANSFER SETTINGS

MiniMed™ 630G Pump to MiniMed™ 770G Pump

## Transfer Settings: MiniMed 630G pump to MiniMed 770G pump

**STEP 1** – Copy your current settings onto the Insulin Pump Settings Form. If for any reason you cannot copy the settings, please contact your healthcare professional.

**STEP 2** – Insert a AA battery into your new pump. You can use either a lithium, alkaline, or rechargeable battery.

**STEP 3** – Follow the instructions found in this guide to enter the settings that you have recorded on the Insulin Pump Settings Form into your new pump. For additional information on your pump, please consult the MiniMed pump system user guide.

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

## Settings from the MiniMed 630G insulin pump

We recommend that you record all settings on this form to refer to as you enter them into the new pump.

### Insulin Settings, Utilities, Sensor Settings

#### INSULIN SETTINGS

BOLUS WIZARD SETUP:

**Menu > Insulin Settings > Bolus Wizard Setup**

Bolus Wizard:  On  Off

Carb Ratio:

12:00 A \_\_\_\_\_ g/U  
 \_\_\_\_\_ g/U  
 \_\_\_\_\_ g/U

Insulin Sensitivity Factor:

12:00 A \_\_\_\_\_ mg/dL/U  
 \_\_\_\_\_ mg/dL/U  
 \_\_\_\_\_ mg/dL/U

BG Target:

12:00 A \_\_\_\_\_ - \_\_\_\_\_ mg/dL  
 \_\_\_\_\_ - \_\_\_\_\_ mg/dL  
 \_\_\_\_\_ - \_\_\_\_\_ mg/dL

Active Insulin Time: \_\_\_\_\_ hr

BASAL PATTERN SETUP:

**Menu > Insulin Settings > Basal Pattern Setup**

Basal 1:

12:00 A \_\_\_\_\_ U/hr  
 \_\_\_\_\_ U/hr  
 \_\_\_\_\_ U/hr

Basal Pattern Name: \_\_\_\_\_

12:00 A \_\_\_\_\_ U/hr  
 \_\_\_\_\_ U/hr  
 \_\_\_\_\_ U/hr

PRESET TEMP SETUP:

**Menu > Insulin Settings > Preset Temp Setup**

Preset Temp Name: \_\_\_\_\_

Type:  Rate \_\_\_\_\_ U/hr  Percent \_\_\_\_\_ %

Duration \_\_\_\_\_ hr

DUAL/SQUARE WAVE:

**Menu > Insulin Settings > Dual/Square Wave**

Dual:  On  Off

Square:  On  Off

#### INSULIN SETTINGS (CONT.)

PRESET BOLUS SETUP:

**Menu > Insulin Settings > Preset Bolus Setup**

Preset Bolus Name: \_\_\_\_\_

Bolus: \_\_\_\_\_ U Type: \_\_\_\_\_

\_\_\_\_\_ Now \_\_\_\_\_ Square

Duration \_\_\_\_\_ hr

BOLUS INCREMENT:

**Menu > Insulin Settings > Bolus Increment**

0.025 U  0.05 U  0.1 U

MAX BASAL/BOLUS:

**Menu > Insulin Settings > Max Basal/Bolus**

Max Basal: \_\_\_\_\_ U/hr Max Bolus: \_\_\_\_\_ U

EASY BOLUS:

**Menu > Insulin Settings > Easy Bolus**

On  Off Step Size \_\_\_\_\_ U

AUTO SUSPEND:

**Menu > Insulin Settings > Auto Suspend**

Alarm:  On  Off

Time: \_\_\_\_\_ hr

BOLUS SPEED:

**Menu > Insulin Settings > Bolus Speed**

Standard  Quick

#### UTILITIES

DISPLAY OPTIONS:

**Menu > Utilities > Display Options**

Brightness \_\_\_\_\_  Auto

Backlight \_\_\_\_\_ sec / min

REMOTE BOLUS:

**Menu > Utilities > Remote Bolus**

On  Off

CARB UNIT:

**Menu > Utilities > Carb Unit**

Grams  Exchanges

#### SENSOR SETTINGS

**Menu > Sensor Settings**

Sensor:

On  Off

### Sensor Settings

HIGH SETTINGS <i>Menu &gt; Sensor Settings &gt; High Settings</i>				
Time Period	High Limit	Alert before High	Alert on High	Rise Alert <sup>a</sup>
12:00 A - _____	_____ mg/dL	<input type="checkbox"/> On _____ min	<input type="checkbox"/> On	<input type="checkbox"/> On _____ mg/dL/min or ↑   ↑↑   ↑↑↑
_____ - _____	_____ mg/dL	<input type="checkbox"/> On _____ min	<input type="checkbox"/> On	<input type="checkbox"/> On _____ mg/dL/min or ↑   ↑↑   ↑↑↑
_____ - _____	_____ mg/dL	<input type="checkbox"/> On _____ min	<input type="checkbox"/> On	<input type="checkbox"/> On _____ mg/dL/min or ↑   ↑↑   ↑↑↑
High Snooze: _____ minutes				

<sup>a</sup> Option only available when Rise Alert is on

LOW SETTINGS <i>Menu &gt; Sensor Settings &gt; Low Settings</i>				
Time Period	Low Limit	Alert before Low	Suspend on Low	Alert on Low
12:00 A - _____	_____ mg/dL	<input type="checkbox"/> On	<input type="checkbox"/> On	<input type="checkbox"/> On
_____ - _____	_____ mg/dL	<input type="checkbox"/> On	<input type="checkbox"/> On	<input type="checkbox"/> On
_____ - _____	_____ mg/dL	<input type="checkbox"/> On	<input type="checkbox"/> On	<input type="checkbox"/> On
Low Snooze: _____ minutes				


## Programming the MiniMed 770G insulin pump

### Startup Wizard

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

Once you've inserted the battery, the pump powers up and the **Startup Wizard** begins, it will guide you through setting language, time format, and current time and date.

1. Select **English**.
2. When the **Select Time Format** screen appears, select a **12 Hour** or a **24 Hour** time format.
3. When **Enter Time** screen appears, the hour will be flashing. Press **^** or **v** to the correct hour and press **o**.
4. The minutes will be flashing. Press **^** or **v** to the correct minutes and press **o**.
5. The **AM / PM** will be flashing. Press **^** or **v** if needed and press **o**. Select **Next**.
6. When the **Enter Date** screen appears, select **Year**. Press **^** or **v** to the correct year and press **o**. Repeat this step to set up **Month** and **Day**. Select **Next**.
7. The rewinding screen will briefly appear, followed by a confirmation message. Select **OK**.




### Setting Max Basal

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

#### Menu > Options > Delivery Settings

1. Select **Max Basal/Bolus** and press **o**.
2. Select **Max Basal** to enter the **U/hr** amount. Use **^** and **v** to set the **Max Basal** amount. Press **o**.
3. Select **Save**.




### Setting Max Bolus

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

#### Menu > Options > Delivery Settings

1. Select **Max Basal/Bolus** and press **o**.
2. Select **Max Bolus** to enter the **U/hr** amount. Use **^** and **v** to set the **Max Bolus** amount. Press **o**.
3. Select **Save**.



### Setting Basal Pattern

**Note:** The pump is intended to be used with a basal pattern. The basal pattern must be manually entered and saved into the pump. There will be no reminder message to program basal rates. Consult a healthcare professional to determine what basal pattern is needed. For more information about basal patterns, consult the MiniMed 770G System User Guide.

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

Menu > Basal > Delivery Settings > Basal Pattern Setup



1. Select **Basal 1**.
2. Select **Options**, then select **Edit**.
3. Press  $\odot$  to enter **End** time. Press  $\odot$  to leave the **End** time at **12:00A**. If you are setting multiple rates use  $\wedge$  or  $\vee$  to change the **End** time and press  $\odot$ .
4. Press  $\wedge$  to enter **U/hr** and press  $\odot$ .
5. Continue adding **End** times and basal rates as needed. The **End** time for your last rate must be **12:00A**.
6. The **Done** option appears only when the last **End** time in your basal pattern is set to **12:00A**. Select **Done** after all settings are entered to go to the next screen.

**Note:** Your basal rates have not yet been saved.

A screen appears that lets you review your basal pattern. Press  $\vee$  to review all the settings. If you need to make any changes, press  $\wedge$  to return to the previous screen (Edit screen).

**Note:** If you press  $\wedge$  to return to the Edit screen without saving, your changes will not be saved. After making any changes on the Edit screen, select **Done** to return to the review screen.

7. Select **Save**. If you do not select **Save**, your changes are not saved.



**CAUTION:** If you have not pressed **Save** after settings are entered and the screen goes dark, the entered settings will not be saved.








**Turn Bolus Wizard On and Setup**

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

**Menu > Options > Delivery Settings > Bolus Estimate Setup**

1. Select **Bolus Wizard** to turn **On**.
2. Press **✓** to see the rest of the message then select **Next**.
3. Review the description of **Carb Ratio** and select **Next**.
4. Press **○** to enter **End** time. Press **○** again to leave the **End** time at **12:00A**. If you are setting multiple **Carb Ratios** use **▲** or **▼** to change the **End** time, then press **○**.
5. Press **▲** or **▼** to enter **g/U** then press **○**. Continue adding **End** times and **Carb Ratios** as needed.
6. Select **Next** when done.
7. Review the description of **Sensitivity Factor** and select **Next**.
8. Press **○** to enter **End** time. Press **○** again to leave the **End** time at **12:00A**. If you are setting multiple **Sensitivity Factors**, use **▲** or **▼** to change the **End** time, then press **○**.
9. Press **▲** or **▼** to enter **mg/dL per U**, then press **○**. Continue adding **End** times and **Sensitivity Factors** as needed.
10. Select **Next** when done.
11. Review the description of **BG Target** and select **Next**.
12. Press **○** to enter **End** time. Press **○** again to leave the **End** time at **12:00A**. If you are setting multiple **BG Targets**, use **▲** or **▼** to change the **End** time, then press **○**.
13. Press **▲** or **▼** to enter **Lo** and then press **○**.
14. Press **▲** or **▼** to enter **Hi** and then press **○**. Continue adding **End** times and **BG Targets** as needed.
15. Select **Next** when done.
16. Review the description of **Active Insulin Time** and select **Next**.
17. Press **○** to enter **Duration** time.
18. Press **▲** or **▼** to enter **hr** then press **○**.
19. Select **Save**.

The **Bolus Wizard** setup is now complete.





**Turn Sensor On**

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

**Menu > Options > Utilities > Sensor Settings**

- Select **Sensor** to turn feature **On**. You can now access the **SmartGuard** menu and enter the settings.





## Programming High Settings

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

### Menu > Options > SmartGuard

1. Select **High Setup**.
2. Press **○** on the time segment. Press **○** again to leave the **End** time at **12:00A**. If setting multiple time segments, press **▲** to enter **End** time of first segment, and press **○**.
3. Press **▲** or **▼** to set **Hi** limit and press **○**. Press **○** to continue.
4. Select each feature you wish to turn on. If a feature is on, select it again to turn it back off.
5. Once settings are selected, select **Next** and then select **Done**.
6. Verify that settings are correct and select **Save**.
7. If snooze time needs to be changed, press **▼** to **Snooze** and press **○**. Select **High Snooze**. Press **▲** or **▼** to the correct time and press **○**. Select **Save**.  
Your **High Setup** is now complete.






## Programming Low Settings

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

### Menu > Options > SmartGuard

1. Select **Low Setup**.
2. Press **○** on the time segment. Press **○** again to leave the **End** time at **12:00A**. If setting multiple time segments, press **▲** to enter **End** time of first segment, and press **○**.
3. Press **▲** or **▼** to set **Lo** limit and press **○**. Press **○** to continue.
4. Select each feature you wish to turn on. If a feature is on, select it again to turn it back off.
5. Once settings are selected, select **Next** and then select **Done**.
6. Verify that settings are correct and select **Save**.
7. If snooze time needs to be changed, press **▼** to **Snooze** and press **○**. Select **Low Snooze**. Press **▲** or **▼** to the correct time and press **○**. Select **Save**.  
Your **Low Setup** is now complete.




Go to *Pair the Transmitter*, page 7 to continue programming sensor settings.

## Turn Auto Mode On

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

### Menu > Options > SmartGuard

1. Select **Auto Mode**.
  2. Select **Auto Mode** again to turn Auto Mode **On**.
  3. Check the screen to make sure that Auto Mode is set to **On**. Select **Save**.
- Note:** If SmartGuard Auto Mode is not yet ready, after selecting **Save** an alert occurs instructing to check the Auto Mode Readiness screen.




**Note:** The SmartGuard Auto Mode warm-up period begins the midnight after the pump starts delivering insulin. If you turn Auto Mode on before the warm up is complete, you may receive an alert at midnight to Enter a BG for Auto Mode. You may want to wait until day 3 of your pump wear to avoid this alert at midnight.


## Pair the Transmitter

**Note:** The settings shown in this guide are for illustration purposes only — your settings may be different.

### Menu > Options > Utilities > Device Options > Pair Device

1. Attach Guardian™ Link (3) transmitter to the charger and place next to the pump.
2. Select **Search** on the pump and immediately remove the transmitter from the charger.
3. The **Select Device** screen appears with a list of available devices. Select the **CGM** device that matches the serial number (**SN**) on your transmitter.
4. Make sure the serial number on the pump screen and transmitter match and select **Confirm**.

Your pump displays a message if the pump and transmitter are paired successfully. If the Sensor feature is turned on, the Connection icon  appears on the Home screen.



**Note:** Please take some time to review ALL your entries and ensure they have been correctly saved.

## Technical Support

### In Need of Technical Assistance?

24-Hour Technical Support is here to help

Call 1.800.646.4633

For medical advice, contact your healthcare team.

# Medtronic



**Medtronic MiniMed**  
18000 Devonshire Street  
Northridge, CA 91325  
USA  
1 800 646 4633  
+1 818 576 5555  
[www.medtronicdiabetes.com](http://www.medtronicdiabetes.com)

RF: M994838A001

© 2022 Medtronic  
M037254C001\_1  
2022-10-19



\*M037254C001\*