English

Indications for Use

The Extended Reservoir is indicated for the subcutaneous infusion of insulin from compatible Medtronic insulin pumps and infusion sets. Refer to your Medtronic insulin pump user guide for a list of compatible insulins and infusion sets.

Contraindications

This reservoir is contraindicated for the infusion of blood or blood products.

Warnings

- Do not use with two or three day infusion sets. The reservoir may only be used for up to seven days with use of the Medtronic Extended infusion set. Use of the reservoir with two or three day infusion sets may lead to hyperglycemia or diabetic ketoacidosis.
- Do not use the reservoir for more than seven days. Using the reservoir for more than seven days may result in the delivery of too little insulin, which may cause hyperglycemia.
- Reservoir and transfer guard are sterile, nonpyrogenic, and for single use only.
- Do not clean or resterilize. Reuse of the reservoir may lead to any one of the following: insulin degradation, infection, inaccurate medication delivery, or leaks which may cause damage to the pump.
- Do not use if package has been opened or damaged.
- Dispose of transfer guard safely in sharps container.

- Failure to remove trapped air from reservoir may result in inaccurate delivery of medication.
- If insulin, or any liquid, gets inside the tubing connector, it can temporarily block the vents that allow the pump to properly fill the infusion set. **This may result in the delivery of too little or too much insulin, which can cause hyperglycemia or hypoglycemia.** If this occurs, start over with a new reservoir and infusion set.
- The pump reservoir should be used for a maximum of seven days, or according to the insulin labeling, whichever duration is shorter.

Precautions

- Read all instructions carefully before using the Extended Reservoir.
- Insulin should be at room temperature before filling reservoir.
- Once reservoir is filled, use promptly. Do not store.
- Check for leaks after changing the reservoir and infusion set.
- Conduct your first setup in the presence of your healthcare professional.

Warranty

Visit <http://www.medtronicdiabetes.com/warranty>.
© 2022 Medtronic. All rights reserved. Medtronic, the Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.

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