

## GLUCOSE MONITORING

### Schedule for adjusting pump settings

When first starting pump therapy or any time pump settings need adjusting:

- Check your blood glucose (BG):
  - When you wake up
  - Before each meal
  - 2 hours after each meal
  - Bedtime
  - Mid-sleep or every 3–4 hours during sleep
- Do not eat between meals.



Checking BGs at these times provides the information needed to adjust and fine-tune pump settings as directed by your healthcare professional.

### Schedule for routine monitoring

Once your pump settings are adjusted correctly and your glucose levels are stable, establish a routine that includes always checking your BG:

- When you wake up
- Before each meal
- Bedtime
- Occasionally mid-sleep
- Test more frequently during travel, times of stress, and illness



## TREATING LOW BLOOD GLUCOSE LEVELS

### How to treat mild/moderate lows

#### 15–15 Rule

If BG drops below 70 mg/dL:

- 1) Eat 15 grams of fast-acting carbohydrate.
- 2) Re-check BG in 15 minutes.
- 3) If BG is still below 70 mg/dL, repeat Steps 1 & 2 every 15 minutes until BG is within range.

#### Items that contain 15 grams:

- 3 to 4 glucose tablets
- 5 jelly beans
- 4 oz juice or soda (not diet)
- 8 oz milk (low or non-fat)
- 1 Tbsp sugar or honey
- If BG is lower than 50 mg/

**If BG is lower than 50 mg/dL, start treatment by eating 20 to 30 grams of carbohydrate or as otherwise directed by your healthcare professional.**

### How to treat a severe low

Keep a Glucagon Emergency Kit on hand in case a severe low occurs. Glucagon can be given by injection to raise glucose levels if you are unable to eat or drink to treat a low, or if you are unconscious.



A family member, co-worker, or friend should be instructed on how to give glucagon.



**NOTE:** If you are using continuous glucose monitoring (CGM), do not rely on sensor glucose values for making treatment decisions or the Suspend on low feature to prevent or treat a low blood glucose.

## TREATING HIGH GLUCOSE LEVELS

### General Guidelines: If BG is high but is lower than 250 mg/dL

- 1) Enter the BG reading into your pump.
- 2) Allow the Bolus Wizard® feature to calculate the correction bolus amount.
- 3) Confirm the bolus amount and select Deliver Bolus.
- 4) Recheck your BG in one hour and again each hour until your BG is back within target range.

**Never ignore high BG readings. Always consult the Bolus Wizard to see if a correction bolus should be taken.**

### General Guidelines: If BG is higher than 250 mg/dL- check for ketones

#### If ketone test is negative:

- 1) Enter BG into pump/consult Bolus Wizard to see if correction dose is needed.
  - Use pump to give correction dose
- 2) Recheck BG in 1 hour
  - If BG is starting to decrease, continue to monitor until normal.
  - If BG is same or higher:
    - Give correction dose using a syringe.
    - Change infusion site, infusion set, reservoir, and insulin.
    - Continue to check BG every hour until BG returns to normal.

#### If ketone test is positive:

- 1) Take correction dose using a syringe.
- 2) Change infusion site, infusion set, reservoir, and insulin.
- 3) Troubleshoot pump.
- 4) Check BG every 1 to 2 hours. Give correction boluses as needed.
- 5) Drink non-carbohydrate fluids.
- 6) If BG continues to rise or if you have moderate to high ketones, nausea, vomiting, or difficulty breathing, notify physician or go to the nearest emergency room.

## DKA PREVENTION

### Sick day Guidelines

**Illness and/or infection usually cause BGs to run higher than normal. Therefore, the risk of developing DKA is increased when you are ill.**

Because DKA symptoms are similar to flu and stomach virus symptoms, check your BG and monitor for ketones often during illness.

- Check BG every 2 hours or as directed by your healthcare professional.
- Check urine or blood for ketones as directed by your healthcare professional.
- Immediately check ketones if you have nausea, vomiting, or abdominal pain.
- Notify doctor if ketones are positive, if you are unable to keep food down, or if no improvement within a few hours. Give a correction dose of insulin with a syringe according to your healthcare professional's recommendations and change infusion set and reservoir.

### Check for ketones

Follow the instructions in your ketone testing kit.



**Unexplained highs that do not decrease with a correction bolus may be caused by a dislodged or kinked infusion set or a weak vial of insulin.**