



**Medtronic**

# Getting Started

with Continuous Glucose Monitoring

## Learning Guide

Paradigm REAL-Time™  
**Revel**  
MINIMED

# Getting Started with Continuous Glucose Monitoring

Keys to Success


# Getting Started with CGM:

## Learning Guide



### TRENDS

CGM finally lets you “fill in the blanks” between your fingersticks.

Focus more on the \_\_\_\_\_ with CGM, instead of each individual number.

CGM tells you about the \_\_\_\_\_ and \_\_\_\_\_ of glucose changes that fingersticks are unable to show.



\_\_\_\_\_ arrow means glucose has gone up or down 1-2 mg/dL per minute.

\_\_\_\_\_ arrows mean glucose has gone up or down over 2.0 mg/dL per minute.



### SG ≠ BG



**BG** - the Blood Glucose readings that are measured by your meter.

**SG** - the Sensor Glucose readings that appear on the CGM display.

You should expect that SG readings on the CGM monitor will be \_\_\_\_\_ than BG meter values.

This is \_\_\_\_\_.

Your meter measures glucose in the \_\_\_\_\_.

Your sensor measures glucose in the \_\_\_\_\_ between the cells (interstitial fluid).

Since glucose is always moving between these two places, there will be a natural difference in BG and SG readings.

# Learning Guide

Most of the time, SG and BG readings will be very \_\_\_\_\_ to each other.

You will see a \_\_\_\_\_ difference between your SG and BG readings when glucose is changing rapidly.

This is likely to occur:


- AFTER a \_\_\_\_\_ or \_\_\_\_\_
- When you see \_\_\_\_\_ or \_\_\_\_\_ arrows on the screen

Always rely on \_\_\_\_\_ readings to make therapy adjustments.

**NOTE:** *Glucose will often change rapidly after treating a low or correcting a high BG.*



## CALIBRATION



You will need to enter BG readings into your pump to ensure your sensor works properly. This is called calibration.

It is best to calibrate \_\_\_\_\_ to \_\_\_\_\_ times a day.

Calibration should be done when blood glucose and interstitial glucose levels are most likely to be \_\_\_\_\_.

This is usually:

- Before \_\_\_\_\_ and \_\_\_\_\_
- When there are no \_\_\_\_\_ or \_\_\_\_\_ arrows on the display.

Calibrating during these times helps keep BG and SG readings close together.


*Calibrate before bed to avoid getting a calibration alarm during the night!*

**ZZZ**

**TIP:** *If BG and SG readings are very different and you don't see arrows on your pump screen, wash your hands, test your BG, and calibrate the system*



## INSERTION



Sensors should be inserted at a \_\_\_\_\_ to \_\_\_\_\_ degree angle.

Inserting the sensor at \_\_\_\_\_ than \_\_\_\_\_ may cause discomfort and/or bleeding.

Inserting the sensor \_\_\_\_\_ than \_\_\_\_\_ degrees can also cause discomfort and may keep the sensor from working properly.

### **REMEMBER: Rotate insertion site with each new sensor**

#### **Choose sites:**

- where your body does not naturally bend a great deal
- where you do not have scars, thickened skin or stretch marks
- that are at least 2 inches away from your infusion set



## WETTING

Sensors need \_\_\_\_\_ to properly wet.

If the sensor is not wetting properly, you may receive

\_\_\_\_\_ and \_\_\_\_\_

alarms during the 2 hours initialization period.

If these alarms occur, *do not worry*. If Sensor Error occurs, simply

clear it. If Lost Sensor occurs within 20 minutes of a sensor start, go to

New Sensor. If the alert occurs after 20 minutes, go to Reconnect Old Sensor.



**If these errors continue to occur during the initialization period, turn the sensor off. After 2 hours have passed, turn sensor on and go to [Link to Sensor and New Sensor](#).**



## PERSONALIZE



High and Low glucose alerts are NOT the same as glucose targets. Sensor \_\_\_\_\_ and \_\_\_\_\_ settings should be set to levels that work best for you.

It is best to start with these alerts set \_\_\_\_\_ to avoid getting too many alerts.

You can then adjust your alert settings over time as you get used to CGM and your glucose control improves.




**If you are getting too many alerts, your settings and repeat times may need to be changed.**



## CARELINK™ SOFTWARE

CareLink software helps you to understand the effects of \_\_\_\_\_, \_\_\_\_\_ and \_\_\_\_\_ on your glucose levels.

CareLink software makes it easier to see areas where changes can be made to improve your glucose control.



### **It is best to start off with these two reports first:**

The \_\_\_\_\_ Report gives you specific information about a particular day including; insulin taken, carbohydrates eaten, sensor readings and alarms. It is helpful to understand how diet, insulin, and activity affect your glucose levels.

The \_\_\_\_\_ Report shows you trends and patterns. These can help you adjust your sensor alert settings and make necessary changes to the amount of insulin you receive.



## LEARNING

A solid understanding of CGM will improve your glucose control and allow you to live your life your way!

However, like anything new, learning CGM takes \_\_\_\_\_ and \_\_\_\_\_.



### **We have many resources to help you along the way.**

- Getting Started with Continuous Glucose Monitoring Guide
- Getting Started with CareLink Software for CGM Guide
- CareLink Software
- 24-Hour HelpLine at: \_\_\_\_\_
- [www.medtronicdiabetes.com](http://www.medtronicdiabetes.com)

**NOTES:** \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

[www.medtronicdiabetes.com](http://www.medtronicdiabetes.com)

**Medtronic Diabetes**  
18000 Devonshire Street  
Northridge, CA 91325-1219  
1.800.646.4633

