

REAL-Time Continuous Glucose Monitoring Values

24 hour Graph
To view your current glucose and the most recent 24 hours of data, press **ACT** twice from the home screen.




3 hour Graph
To view your current glucose and the most recent three hours of data, press **ACT** once from the home screen.



Pump Modes

Normal



The pump is ON and delivering as programmed with no special circumstances.

Special



The pump is ON and delivering as programmed, but there is a special circumstance, such as Dual/Power OFF, Prime, Alarms, or Suspended, Empty Reservoir, or one of the following reasons:
The pump is not delivering for

Attention



Insulin Delivery Stopped.
The pump is communicating with the transmitter.

Sensor Communicating



The Sensor Feature is turned ON, but the pump is not communicating with the transmitter.

Sensor On



THE MINIMED
Paradigm REAL-Time
INSULIN PUMP AND CONTINUOUS GLUCOSE MONITORING SYSTEM

Menu Map




Medtronic MiniMed
Northridge, CA 91325
USA
800-646-4633 (800-MiniMed)
818-576-5555

6025201-012

111406



Menus and Menu Options



Scroll bar (appears when more screen text is available)

Press to scroll down to view additional screen text.

Press to scroll back up.

Press the ACT button from the Main Menu to get to the sub-menus



BOLUS MENU

- Set Bolus¹
 - ▶ Set Bolus
- Bolus History
- Bolus Wizard Setup
 - ▶ Edit Settings
 - ▶ Wizard
 - ▶ Carb Units
 - ▶ Carb Ratios
 - ▶ BG Units
 - ▶ Sensitivity
 - ▶ BG Target
 - ▶ Active Ins Time
 - ▶ Review Settings
- Max Bolus
 - ▶ Max Bolus Setup
- Dual/Square Bolus
 - ▶ Off
 - ▶ On
- Easy Bolus
 - ▶ Off
 - ▶ On/ Set
- BG Reminder
 - ▶ Off
 - ▶ On

1. Displays only when the Bolus Wizard[®] feature is off.

SUSPEND

Press ACT to Stop Pump

SENSOR MENU

- Enter Meter BG²
- Sensor Update Hist²
- Sensor Alarm History²
- Sensor Start²
 - ▶ New Sensor
 - ▶ Reconnect Old Sensor
 - ▶ Find Lost Sensor
- Sensor Setup
 - ▶ Edit Setting
 - ▶ Review Settings

2. Displays if Sensor is ON.

BASAL MENU

- Set/Edit Temp Basal
 - ▶ Set Duration
- Cancel Temp Basal³
- Select Patterns⁴
- Set/Edit Basal
 - ▶ Set Basal Rate 1
- Basal Review
- Max Basal Rate
- Patterns
 - ▶ Off
 - ▶ On
- Temp Basal Type
 - ▶ Set Temp Basal As
 - ▶ Insulin Rate
 - ▶ Percent of Basal

3. Displays only after setting Temp Basal.

4. Displays only when the Patterns feature is ON.

PRIME MENU

- Manual Prime⁵
- Fixed Prime
- Rewind
- Prime History

5. Displays only after a rewind.

UTILITIES MENU

- Lock Keypad
- Alarm
 - ▶ Alarm History
 - ▶ Alert Type
 - ▶ Auto-Off
 - ▶ Low Resv Warning
 - ▶ Insulin Units
 - ▶ Time
- Daily Totals
- Time / Date
 - ▶ 12 Hour Setup
 - ▶ 24 Hour Setup
- Alarm Clock
 - ▶ Off
 - ▶ On/ Set
- Meter Options
 - ▶ Off
 - ▶ On
- Remote Options
 - ▶ Off
 - ▶ On
- Block
 - ▶ Off
 - ▶ On
- Selftest
- User Settings
 - ▶ Save Setting
 - ▶ Restore Settings
 - ▶ Clear Settings
 - ▶ History
- Language



Press the ESC button from from a submenu screen to return to the Main Menu

Note: All screens are samples only. Actual screens depend on the current active functions. Blue Arrows (▶) indicate submenus.

Pump Status Screen (Sample Screen Only)

STATUS U100
Last Bolus N 3.8U
9:07A 18 MAR
Basal 1: 0.15U/H Reservoir
Started:
13MAR, 12:44P
Units left: 144.0U
Time left: > 24 hours
BG Reminder in 2:04h
Meter: On
Battery: Normal
Auto Off - 12HR
Fri 21 MAR 2004
S/N# 123456
Paradigm 522
VER X.XX X.X

Sensor Status Screen (Sample Screen Only)

SENSOR STATUS
Next Cal: 08:30
Sensor Age: 2d 3h
Sensor Isig: 25.97
Transmtr Batt: Good
Transmtr Ver: 1.1
Transmtr ID: 2111111
Pump S/N#: 2222222



The ACT button opens menus and activates menu selections.



From the Home screen, the Escape button opens the Status screen or Sensor graph screens. From all other screens, the Escape button moves backward to a previous menu, or exits a screen without setting a value.



The Backlight / Down Arrow button turns on the backlight from the Home screen, scrolls down through the menus, and changes entry values.



From the Home screen, the Easy Bolus button is a shortcut to set and deliver an Easy Bolus. From all other screens, the Up Arrow button scrolls up through the menu, and changes entry values.



The Express Bolus button gives you immediate access to the Bolus or Bolus Wizard feature menus.