

## Part 1. Inserting a New Sensor

Wash your hands and clean insertion site with alcohol.

### 1. Open sensor package.

Pull corner of paper covering to open sensor package.



### 2a. Hold sensor by plastic pedestal.

Remove sensor with attached pedestal by holding pedestal only. Place sensor/pedestal on a clean, flat surface (such as a table).



### 2b. Tuck adhesive tab.

Make sure that sensor's adhesive tab is tucked under the snaps and sensor connector.



Correct



Incorrect

### 3. Load sensor into serter.

Grip serter exactly as shown **with thumb on thumb print on serter. Do not hold green buttons.** Push serter down onto pedestal until base of serter sits flat on table.



Correct



Incorrect




### 4. Detach serter from pedestal.

To detach serter from pedestal, grip serter as shown, with thumb on thumb print on serter. With the other hand, place two fingers on pedestal feet and slowly pull serter straight up.

**NOTE: Make sure that pedestal is firmly on the table before pulling serter away.**

**Warning: Do not detach pedestal from serter in mid-air as this may damage sensor.**

Fingers NOT holding green buttons



### 5a. Place serter on body.

Place base of serter flat against cleaned insertion site. Hold serter steadily against your cleaned insertion site, without pushing serter too deeply into skin.



**NOTE: Failing to hold serter securely flat against the body may allow serter to spring back after pressing the buttons and result in improper insertion of the sensor.**

### 5b. Insert sensor.

Press and release the **bump** on both buttons at same time.



### 5c. Hold serter against body.

Continue holding serter against your body to allow adhesive time to stick to skin.




### 5d. Remove serter from body.

Slowly pull serter away from skin, making sure buttons are not pressed.



### 6. Remove needle housing.

Gently hold base of sensor against skin with one hand. With the other hand, hold needle housing **at the top** and slowly pull straight out, away from sensor. Dispose needle housing in sharps container.



### 7a. Remove adhesive pad liner.

Hold sensor in place and gently remove liner from under adhesive pad.



### 7b. Press entire adhesive pad to skin.

Firmly press adhesive against skin and smooth entire adhesive pad so it sticks to skin.

**NOTE: Enlite adhesive is pressure-sensitive. Continue applying pressure to ensure sensor remains inserted in skin for whole 6 days of wear.**



### 8a. Untuck adhesive tab.

Untuck adhesive tab from under sensor connector.



### 8b. Straighten adhesive tab.

Straighten adhesive tab so it lies flat against your skin, but do not remove adhesive liner yet.




## Part 2. Taping the Sensor

1. Remove liner marked 1 from overtape. Do not remove two smaller liners on sides of overtape.



2. **Important: Attach overtape to both sensor and skin next to sensor.**

Tape is over both sensor and skin



3. Stretch remaining part of overtape around sensor connector so that overtape sticks to curved adhesive pad and does not block sensor connector and snaps. **Continue to press overtape to your skin to help ensure that it sticks securely.**




Sensor connector and snaps in hole of tape

4. Remove two liners marked 2 from the sides of the overtape and press adhesive against the skin.




5. This image is an example of overtape applied correctly. Sensor connector and snaps are not covered and appear in opening of overtape.



## Part 3. Connecting the Transmitter


1. With one hand, hold sensor in place. With the other hand, connect transmitter to sensor.




2. You will hear a faint "click" that indicates the two components are connected. Check for a green light to flash on the transmitter.



3. Remove liner on adhesive tab.



4. Fold adhesive tab over and onto transmitter. **Important: Be careful not to pull adhesive tab too tightly or it may cause transmitter to pull from sensor connector.**





5. Press adhesive onto transmitter.



## Part 4. Initializing the Sensor



1. From Home Screen, press **ACT** to open **Main Menu**.
2. Scroll to **Sensor** and press **ACT**.
3. Scroll to **Link to Sensor** and press **ACT**.
4. Press **ACT** to select **New Sensor**
5. You will see a message that your sensor will be ready in 2 hours.



## Part 5. Calibrating


Two hours after sensor has started, "METER BG NOW" alert will notify you for first calibration. There are two ways to calibrate:

1. Press **CB** button
2. Enter BG meter value
3. You will receive message "BG to update sensor?" Select "Yes" and press **ACT**.



If you are using a linked meter, your BG meter value will automatically appear on your pump screen. Then follow step 3 above.

1. From Home screen, press **ACT** to open **Main Menu**.
2. Scroll to **Sensor Menu** and press **ACT**.
3. Press **ACT** to select **Calibrate** (You will only see the option to Calibrate if the sensor and transmitter are connected and sensor is ready for first calibration).



**NOTE:** Calibrating through Sensor Menu can only occur if pump's insulin delivery is not suspended.