

I am wearing a required medical device prescribed by my physician.



Patient Information

my name and address:

physician name and telephone:

Medical Device Information

medical device type:

device serial number:

This device can withstand exposure to common electrostatic (ESD) and electromagnetic interference (EMI), including airport security systems. During exposure, it can remain connected and powered on.

While onboard aircraft, turn off (power off) any wireless transmitter components or remote controls for this device.

Do not expose this device to any radiation (such as X-ray, MRI, CT scan).



Medtronic

MINIMED

Medtronic MiniMed
Northridge, CA 91325
USA
800-646-4633 (800-MiniMed)
818-576-5555

E.U. Representative
Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
31 (0) 45 566 8000
www.minimed.com

R_x Only

6025051-0U1 060106