

"The ease and simplicity of sharing data with MiniMed® Connect is a game changer, and the data is extremely powerful... MiniMed Connect will absolutely enhance our capability to improve patient care."†

DR. LAURA GANDRUD
Endocrinologist



**BETTER DATA.
BETTER CARE.**

Save time in office visits with automatic uploads to CareLink®

MiniMed Connect links a MiniMed pump and CGM system to an iPhone® or iPod touch®. Your patients' data is automatically sent to CareLink® Personal. Your office can easily access the data – **eliminating the need for manual pump downloads during visits.**

Diabetes data in easy reach for your patients

MiniMed Connect lets your patients check their sensor glucose levels discreetly and conveniently. The more they see their levels, the more they can adjust for better control. Your patients can also allow loved ones to see how they're doing and receive text messages in case they need help.



GET CONNECTED

Talk with your Medtronic representative
or visit medtronicdiabetes.com

† All statements made by featured healthcare professionals reflect their independent professional judgment. Participant was not compensated for the use of their quote and statement. Participant has been a paid consultant in the past. Results may vary.



1

The uploader sends pump and sensor glucose data to the app on an iPhone® or iPod touch® about every five minutes.

2

The app gives patients an easy way to check details, such as:

- Sensor glucose data
- Active insulin level
- Calibration time
- Insulin in reservoir
- Sensor life
- Battery level
- Connection status

3

Data is automatically uploaded to CareLink Personal every 24 hours. Patients can set up automatic text alerts and give friends and family access to information online.

MiniMed Connect works with the MiniMed® 530G and MiniMed® Paradigm® Revel™ systems. No prescription is needed.



The **MiniMed Connect** app for iPhone® and iPod touch® can be downloaded from the App StoreSM. We're currently working on a version for Android™.



Important Safety Information: CareLink®

MiniMed Connect is intended to provide a secondary display of continuous glucose monitoring and/or insulin pump data on a suitable consumer electronic device to care partners and users of a MiniMed 530G system or MiniMed® Paradigm® Revel™ system for the purposes of passive monitoring. MiniMed Connect system is not intended to replace the real-time display of continuous glucose monitor and/or insulin pump data on the primary display device (i.e. the sensor-augmented pump). All therapy decisions should be based on blood glucose measurements obtained from a blood glucose meter.

CareLink Personal therapy management software is intended to be used together with advice from a healthcare professional familiar with the diagnosis and treatment of diabetes. Do not make any changes to treatment without talking to a healthcare professional first.

Please visit medtronicdiabetes.com/important-safety-information for complete details.

MiniMed and CareLink are registered trademarks of Medtronic MiniMed, Inc. Apple, iPhone and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android is a trademark of Google, Inc. Wireless carrier; text and/or data charges may apply. Mobile device not included.

930M14374-012 20160125 © 2016 Medtronic MiniMed, Inc. All rights reserved.

GET CONNECTED

Talk with your Medtronic representative or visit medtronicdiabetes.com

Medtronic