

We recommend that you record all settings on this form to refer to as you enter them into new pump.

INSULIN SETTINGS

BOLUS WIZARD SETUP

Menu > Insulin Settings > Bolus Wizard Setup

Bolus Wizard: On Off

Carb Ratio:

12:00 A _____ - _____ g/U
 _____ - _____ g/U
 _____ - _____ g/U

Insulin Sensitivity Factor:

12:00 A _____ - _____ mg/dL per U
 _____ - _____ mg/dL per U
 _____ - _____ mg/dL per U

BG Target:

12:00 A _____ - _____ mg/dL
 _____ - _____ mg/dL
 _____ - _____ mg/dL

Active Insulin Time: _____ hr

BASAL PATTERN SETUP

Menu > Insulin Settings > Basal Pattern Setup

Basal 1: _____ Basal Pattern Name: _____
 12:00 A _____ U/hr 12:00 A _____ U/hr
 _____ U/hr _____ U/hr
 _____ U/hr _____ U/hr

PRESET TEMP SETUP

Menu > Insulin Settings > Preset Temp Setup

Preset Temp Name: _____
 Type _____ Percent _____ % Duration _____ hr

DUAL/SQUARE WAVE:

Dual: On Off Square: On Off

PRESET BOLUS SETUP

Menu > Insulin Settings > Preset Bolus Setup

Preset Bolus Name: _____
 Bolus: _____ U Type: _____
 % _____ Now _____ Square Duration _____ hr

INSULIN SETTINGS (continued)

PRESET BOLUS SETUP

Menu > Insulin Settings > Preset Bolus Setup

Preset Bolus Name: _____
 Bolus: _____ U Type: _____
 % _____ Now _____ Square Duration _____ hr

BOLUS INCREMENT: 0.025U 0.05U 0.1U

MAX BASAL/BOLUS:

Max Basal: _____ U/hr Max Bolus: _____ U

EASY BOLUS: On Off Step Size _____ U

AUTO SUSPEND: On _____ hr Off

BOLUS SPEED: Standard Quick

REMINDERS

Menu > Reminders

PERSONAL:

Name: _____ Time: _____ AM / PM
 Name: _____ Time: _____ AM / PM

BOLUS BG CHECK: On Off

MISSED MEAL BOLUS:

Meal 1: _____ AM / PM - _____ AM / PM
 Meal 2: _____ AM / PM - _____ AM / PM
 Meal 3: _____ AM / PM - _____ AM / PM

LOW RESERVOIR:

Units: _____ U Time: _____ hr

SET CHANGE:

On _____ days Off

UTILITIES

Menu > Utilities

DISPLAY OPTIONS:

Brightness _____ Auto Backlight _____ sec / min

REMOTE BOLUS: On Off

CARB UNIT: Grams Exchanges

SENSOR SETTINGS: LOW SETTINGS Menu > Sensor Settings > Low Settings

Time Period	Low Limit	Alert before Low	Suspend on Low	Alert on Low
12:00 A - _____	_____ mg/dL	<input type="checkbox"/> On	<input type="checkbox"/> On	<input type="checkbox"/> On
_____ - _____	_____ mg/dL	<input type="checkbox"/> On	<input type="checkbox"/> On	<input type="checkbox"/> On
_____ - _____	_____ mg/dL	<input type="checkbox"/> On	<input type="checkbox"/> On	<input type="checkbox"/> On

Low Snooze: _____ minutes

SENSOR SETTINGS: HIGH SETTINGS Menu > Sensor Settings > High Settings

Time Period	High Limit	Alert before High	Alert on High	Rise Alert
12:00 A - _____	_____ mg/dL	<input type="checkbox"/> On _____ min	<input type="checkbox"/> On	<input type="checkbox"/> On _____ mg/dL/min or ↑ ↑↑ ↑↑↑
_____ - _____	_____ mg/dL	<input type="checkbox"/> On _____ min	<input type="checkbox"/> On	<input type="checkbox"/> On _____ mg/dL/min or ↑ ↑↑ ↑↑↑
_____ - _____	_____ mg/dL	<input type="checkbox"/> On _____ min	<input type="checkbox"/> On	<input type="checkbox"/> On _____ mg/dL/min or ↑ ↑↑ ↑↑↑

High Snooze: _____ minutes