


Description


Name:	i-port Advance®
Detail:	i-port Advance [®] is the first injection port to include an inserter. This simple feature provides the user with a virtually painless insertion when applying the soft cannula (a small, flexible tube) under the skin. The result is a comfortable, dependable application experience. [*] Once applied, only the soft cannula remains under the skin, allowing for direct medication delivery from the injection into the subcutaneous tissue.

Dimensions and Characteristics

Height	0.366 inch (9.3 mm) above skin, applied	Cannula Length	Model Number	NDC
Diameter - body	1.128 inches (28.7 mm)	6 mm	MMT-100	76300-0100-10
Diameter - adhesive pad	1.555 inches (39.5 mm)		10100	/0500 0100 10
Length - soft cannula	Available in 6 mm and 9 mm lengths	9 mm	MMT-101	76300-0101-10
Caliber - insertion needle 27 gauge				
Residual space - delivery channel	< 0.0026 mL (0.26 IU)			
U.S. Availability	By prescription only			
Indications for use	i-port Advance® is indicated for both adults and children requiring daily subcutaneous injections.			
Needle Compatibility	Injection needle should not exceed 5/16" (8mm) in length			
Septum durability	Septum durability Up to 75 needle sticks			
Wear time	Up to 72 hours (or 75 needles sticks, whatever comes first)			
Packaging	Each i-port Advance® comes packaged individually in a box of 10 devices or roughly a 1month supply			



Cross-section: i-port Advance[®] in Use



*Survey data on file; individual results may vary.

Important Safety Information: The i-port Advance injection port is prescription-indicated for patients who administer, or to whom is administered, multiple daily subcutaneous injections of physician prescribed medications, including insulin. The device may remain in place for up to 72 hours to accommodate multiple drug injections without the discomfort of additional needle sticks. Site infection and/or irritation are possible. Patients experiencing such symptoms should immediately consult their healthcare provider and/or discontinue use.

For more information about i-port Advance, please contact:

Medtronic Diabetes at 1-800-646-4633.

www.i-port.com

Medtronic Diabetes 18000 Devonshire Street Northridge, CA 91325

